

Command	ACLEAR
PURPOSE	Clear analyzer spectrum
PARAMETERS	
ANLID	Analyzer identifier, may contain the "*" wildcard to select a group of analyzers.
/LIST	List names of analyzers.
FUNCTION	All spectrum elements are set to zero.
REMARK	Analyzers, for which the property "PROTECTED" has been attributed by the AGEN macro, the AGEN command or the AMOD command, are not affected by this command. This is a useful option for analyzers which store calibration information that is not to be modified during accumulation of different input files.
EXAMPLE	<p>ACLEAR DELTA Clear analyzer DELTA.</p> <p>ACLEAR * Clear all analyzers.</p> <p>ACLEAR MASS(*) Clear all elements of the indicated array of analyzers.</p>